



Contents lists available at ScienceDirect

International Dairy Journal

journal homepage: www.elsevier.com/locate/idaairyj

Short communication

Colorimetry of dairy products

Rehana P. Ponnal^a, Jackie E. Wood^{b,*}, Brendon D. Gill^b, Carlos A. Bergonia^b,
Wendy M. Longstaff^b, Valerie Slabbert^b, Lissa C. Bainbridge-Smith^a, Robert A. Crawford^c^a Centre for Applied Science & Primary Industries, Wintec, Private Bag 3036, Waikato Mail Centre, Hamilton 3240, New Zealand^b Fonterra Waitoa, PO Box 7, Waitoa, New Zealand^c Fonterra Te Rapa, PO Box 10397, Hamilton 3241, New Zealand

ARTICLE INFO

Article history:

Received 24 June 2020

Received in revised form

14 September 2020

Accepted 20 September 2020

Available online 1 October 2020

ABSTRACT

Reflectance colorimetry was investigated to measure the colour of whole (WMP) and skim milk (SKMP) powders, infant formulas (IF), an adult nutritional product (ANP), unsalted (USB) and salted (SB) butters. The colorimeter gave precise results and was cost effective, illustrating the benefits of the adoption of a colorimeter rather than the current practice of using reference charts and powders. New Zealand dairy products are yellower compared with other dairy products because of their relatively high naturally occurring β -carotene content.

© 2020 Elsevier Ltd. All rights reserved.

1. Introduction

Food colour is measured by colour chart comparison, visible near infrared spectroscopy, computer-based image analysis, reflectance spectroscopy, transmittance spectroscopy and colorimetry (Dufossé & Galaup, 2010; Pathare, Opara, & Al-Said, 2013). Advantages of colorimetry are (i) it is performed in uniform light, (ii) it is not reliant on human subjective measurements, (iii) it is portable and cheap, (iv) samples can be solid or an opaque liquid, and (v) it is not reliant on storage of references that may vary with time.

The purpose of this study was to examine whether colorimetry could replace colour charts and reference powder comparisons and to apply this technique to the routine sensory evaluation of dairy products, including adult nutritional powders and stage three infant formulas not previously reported. For adoption to be successful, the colorimeter needs to be easy to use and cost effective and give precise results comparable to the current methods.

2. Materials and methods

The colour of skim milk, whole milk, infant formulas, adult nutritional powders, winter and early spring manufactured salted and unsalted butter were measured using a Konica Minolta CR400

colorimeter (ThermoFisher, Auckland, New Zealand). Daily calibration used the reference white ceramic tile. A Konica Minolta CR A50 (50 × 9 mm glass-topped colorimetry cell) was used for powders and this was packed with 10 g of the powder. The surface was smoothed using a ruler. Between powder sample measurements, the cell was emptied and any residual was removed using brushes.

A fresh flat surface on refrigerated butter was cut with a knife and before being covered with a layer of transparent polyethylene film. To calibrate for butter the reference tile was covered with a single layer of the film.

To assess the repeatability and intermediate precision of the colorimeter, a single batch of adult nutritional powder (stored in gas-flushed and re-sealable foil bags) was used. For repeatability, the powder was measured in duplicate 15 times, while for intermediate precision it was measured in duplicate on 23 days over 2.5 months.

Colorimetry ranges for 60 different batches of skim, whole and different nutritional milk powders, salted and unsalted butters were measured. Butter colour reference tiles were also measured.

The reference powder method is set at point two, using a five-point scale, where good powder is within one point of the reference. Acceptable butter must be the same or be less coloured than the US Department of Agriculture (USDA) medium light reference.

3. Results and discussion

Results are reported in Lab colour space where L^* is lightness from black (0) to white (100), a^* is green (–) to red (+) and b^* is

* Corresponding author. Tel.: +64 7 858 1462.

E-mail address: jackie.wood@fonterra.com (J.E. Wood).

Table 1
Results for repeatability and intermediate precision.^a

Variable	Mean single day	RSD _r (%)	R	Mean 2.5 months	RSD _{ir} (%)
L*	91.08	0.16	0.40	91.00	0.35
a*	-6.85	-1.39	0.26	-6.68	-1.51
b*	17.55	1.73	0.84	17.38	1.34

^a For mean single day, n = 15 duplicate pairs; for mean 2.5 months, n = 23 duplicate pairs. Abbreviations: RSD_r, repeatability relative standard deviation of n pairs of duplicates; R, repeatability limit (= standard deviation × t_∞ × √2); RSD_{ir}, intermediate repeatability relative standard deviation of n pairs of duplicates.

blue (-) to yellow (+) (ISO/CIE, 2008). The results demonstrate good repeatability and intermediate precision of the colorimeter (Table 1). The range and precision for different products are shown in Table 2.

Green colour in milk is due to the riboflavin content (Nozière et al., 2006). Processing can form Maillard browning reaction products increasing redness and this, plus changes in particle size

Table 2
Results for the different product types.^a

Sample	Number/pairs	L*				a*				b*						
		Mean	Range	RSD _r (%)	R	RSD _{ir} (%)	Mean	Range	RSD _r (%)	R	RSD _{ir} (%)	Mean	Range	RSD _r (%)	R	RSD _{ir} (%)
WMP	124/62	89.48	87.41–91.66	0.39	0.96	1.05	-7.94	-6.73 to -8.74	-0.50	0.11	-6.87	22.00	19.21–23.90	0.57	0.34	6.39
Ref WMP	2/0	90.91	91.52–90.29				-6.90	-7.11 to -6.96				22.69	22.95–22.43			
SKMP	136/68	91.44	89.35–92.77	0.39	1.00	0.95	-7.44	-6.88 to -9.30	-1.02	0.21	-7.75	16.63	14.43–21.55	1.24	0.57	9.81
Ref SKMP	12/6	91.81	90.38–92.69				-7.16	-6.91 to -7.52				15.76	15.05–17.30			
IFS1	164/82	87.45	84.45–89.63	0.32	0.77	1.62	-4.39	-0.44 to -5.21	-18.12	2.20	-10.29	18.00	14.44–21.27	0.97	0.48	14.17
Ref IFS1	14/7	87.24	85.98–89.22				-4.39	-3.92 to -4.68				18.88	15.12–20.56			
IFS2	120/60	86.90	85.21–88.83	0.40	0.90	1.07	-4.88	-0.45 to -6.32	-7.48	1.01	-17.26	19.20	16.18–21.21	0.83	0.44	8.26
Ref IFS2	14/7	86.91	85.23–88.51				-5.00	-4.27 to -6.00				19.24	17.22–21.21			
IFS3	120/60	88.36	86.75–89.87	0.36	0.87	0.82	-6.56	-5.11 to -7.77	-0.67	0.12	-15.79	20.20	16.48–23.83	0.64	0.36	9.60
Ref IFS3	18/9	88.86	87.60–89.87				-5.93	-5.30 to -7.42				18.52	16.66–20.23			
ANP	120/60	90.64	87.99–91.61	0.36	0.90	0.64	-6.87	-5.53 to -16.59	-13.35	2.54	-13.62	17.09	15.91–18.00	1.11	0.52	2.23
Ref ANP	14/7	90.52	87.99–91.50				-6.49	-5.53 to -7.03				16.88	16.33–17.63			
USB	120/60	87.82	82.93–90.67	0.70	1.70	1.52	-4.00	-3.23 to -4.73	-4.03	0.45	-6.93	30.37	26.81–33.41	2.04	1.71	4.65
SB	120/60	85.76	81.89–89.77	0.88	2.10	1.28	-3.99	-2.57 to -5.51	-5.25	0.58	-7.39	32.07	25.53–35.06	2.58	2.29	3.90

^a Mean and range of n replicates shown. Abbreviations: RSD_r, repeatability relative standard deviation of n pairs of duplicates; R, repeatability limit (= standard deviation × t_∞ × √2); RSD_{ir}, intermediate repeatability relative standard deviation of n pairs of duplicates.

Table 3
Literature colorimetric results for milk powders, infant formula and butters.^a

Product	Location	L*	a*	b*	References
SMP	USA	79.05	0.78	7.36	Koc, Heinemann, & Ziegler (2003)
SMP	Austria	94.9	-1.7	11.3	Kneifel, Ulberth, & Schaffer (1992)
SMP	New Zealand	93.0	-2.2	16.4	
SMP	USA	92.5	-2.6	18.3	
WMP	Europe	95.6	-3.6	19.8	
WMP		95.96	-3.05	9.56	Nanua, McGregor, & Godbert (2000)
SMP		96.94 ± 0.56	-2.32 ± 0.14	11.12 ± 1.32	Pugliese et al. (2017)
WMP		96.01 ± 0.54	-1.74 ± 0.31	14.45 ± 0.55	
WMP, autumn [†]	Argentina	92.63 ± 0.44	-0.48 ± 0.11	19.61 ± 0.57	Grigorini et al. (2007)
WMP, autumn [‡]		93.38 ± 0.08	-0.73 ± 0.10	18.62 ± 0.42	
WMP, winter [†]		93.09 ± 0.17	-0.53 ± 0.13	20.30 ± 0.66	
WMP, winter [‡]		93.19 ± 0.30	-0.65 ± 0.17	18.55 ± 0.54	
WMP, spring [†]		92.56 ± 0.35	-0.65 ± 0.20	21.11 ± 0.61	
WMP, spring [‡]		93.22 ± 0.18	-0.50 ± 0.19	19.12 ± 0.66	
WMP, summer [†]		92.19 ± 0.32	-0.74 ± 0.10	20.37 ± 0.58	
WMP, summer [‡]		92.90 ± 0.48	-0.85 ± 0.16	20.05 ± 0.92	
WMP	United Arab Emirates	95.37	-0.01	12.67	Suliman, Elamin, Elkhalfifa, & Laleye (2014)
Infant formula, stage 1	Europe	92.34	-0.25	12.31	Ferrer, Alegria, Farré, Clemente, & Calvo (2005)
Infant formula, stage 1		90.70	-0.25	14.78	
Butter 10 °C, summer		63.2	3.9	31.1	Kneifel et al. (1992)
Butter 10 °C, winter		70.8	4.7	28.3	

^a Abbreviations are: SMP, skim milk powder; WMP, whole milk powder. Heat treatments indicated as: [†], indirect; [‡], direct.

affecting light diffusion, changes lightness (Grigioni et al., 2007; Pugliese et al., 2017). Previous literature measurements of L*a*b* scores are given in Table 3. Powders in this study were darker than powders previously reported.

New Zealand whole milk powders and butters have high yellowness scores from β-carotene in grass feed (Dufossé, Fernández-López, Galaup, & Pérez-Alvarez, 2015). Skim milk powder is usually darker than whole milk because of light diffusion, but because of the yellowness whole milk powder was darker (Pugliese et al., 2017).

Infant formulas were generally darker and less green than the milk powders or European infant formulas. Change in colour was probably due to the extra minerals and oils ingredients used to manufacture infant formula. All butters were within the light to medium range compared to USDA butter colour reference tiles reported in Table 4. The colorimeter was more cost effective based on the 5 min allowed for chart and powder comparisons whereas colorimetry measurements took less than 1 min.

Table 4
USDA butter colour guide measurements.

Plate colour description	<i>L</i> *	<i>a</i> *	<i>b</i> *
Very light (VL)	92.03	-2.44	22.43
Light (L)	90.91	-2.39	25.14
Medium light (ML)	89.19	-2.26	31.30
Medium (M)	87.86	-1.38	36.27
Medium high (MH)	86.53	0.27	37.19
High (H)	84.97	2.67	38.95
Very high (VH)	77.54	9.10	44.22

4. Conclusions

Colorimetry is a feasible alternative to using reference colour charts or powder references as it gives precise results and is more cost effective. New Zealand dairy products differ in colour from similar international products previously reported, because of their significant yellow colour from the high naturally occurring β -carotene content.

CRediT author statement

Rehana P. Ponnal: Validation, formal analysis, investigation, data curation.

Jackie E. Wood: Conceptualization, methodology, validation, formal analysis, data curation, writing original draft, writing and editing, supervision, resources, project administration, visualisation.

Brendon D. Gill: Software, writing original draft, writing and editing, visualisation.

Carlos A. Bergonia: Conceptualization, methodology, resources, project administration, supervision.

Wendy M. Longstaff: Conceptualization, methodology, resources, supervision.

Valerie Slabbert: Resources, project administration.

Lissa C. Bainbridge-Smith: Supervision and project administration.

Robert. A. Crawford: Conceptualization, methodology, writing and editing, supervision, project administration.

Acknowledgements

The authors would like to acknowledge Alan Russell and Marie Edlington of ThermoFisher for the loan of the CR400 colorimeter for this project and Steve Holroyd and Claire Woodhall for proof-reading this manuscript.

References

- Dufossé, L., Fernández-López, J., Galaup, P., & Pérez-Alvarez, J. A. (2015). Color measurements of muscle-based and dairy foods. In L. M. L. Nollet, & F. Toldrá (Eds.), *Handbook of food analysis* (3rd ed., Vol. 1, p. 3–20). Boca Raton, FL, USA: CRC Press.
- Dufossé, L., & Galaup, P. (2010). Color of dairy foods. In L. M. L. Nollet, & F. Toldrá (Eds.), *Handbook of dairy foods analysis, Part IV: Sensory quality* (pp. 581–602). Boca Raton, FL, USA: CRC Press.
- Ferrer, E., Alegría, A., Farré, R., Clemente, G., & Calvo, C. (2005). Fluorescence, browning index, and color in infant formulas during storage. *Journal of Agricultural and Food Chemistry*, 53, 4911–4917.
- Grigioni, G., Biolatto, A., Irurueta, M., Sancho, A. M., Páez, R., & Pensel, N. (2007). Colour changes of milk powder due to heat treatments and seasons of manufacture. *Ciencia y Tecnología Alimentaria*, 5, 335–339.
- ISO/CIE. (2008). *ISO 11664-4:2008(en)/CIE S 014-4/E:2007. Joint ISO/CIE standard: Colorimetry—Part 4: CIE 1976 L*a*b* colour space*. Geneva, Switzerland: International Organisation for Standardisation.
- Kneifel, W., Ulberth, F., & Schaffer, E. (1992). Tristimulus colour reflectance measurement of milk and dairy products. *Le Lait*, 72, 383–391.
- Koc, A. B., Heinemann, P. H., & Ziegler, G. R. (2003). A process for increasing the free fat content of spray-dried whole milk powder. *Journal of Food Science*, 66, 210–216.
- Nanua, J. N., McGregor, J. U., & Godbert, J. S. (2000). Influence of high-oryzanol rice bran oil on the oxidative stability of whole milk powder. *Journal of Dairy Science*, 83, 2426–2431.
- Nozière, P., Graulet, B., Lucas, A. E., Martín, B. C. R., Grolier, P., & Doreau, M. T. (2006). Carotenoids for ruminants: From forages to dairy products. *Animal Feed Science and Technology*, 131, 418–450.
- Pathare, P. B., Opara, U. L., & Al-Said, F. A.-J. (2013). Colour measurement and analysis in fresh and processed foods: A review. *Food and Bioprocess Technology*, 6, 36–60.
- Pugliese, A., Cabassi, G., Chiavaro, E., Paciulli, M., Carini, E., & Mucchetti, G. (2017). Physical characterization of whole and skim dried milk powders. *Journal of Food Science & Technology*, 54, 3433–3442.
- Suliman, A. M. E., Elamin, O. M., Elkhalfia, E. A., & Laleye, L. (2014). Comparison of physicochemical properties of spray-dried camel's milk and cow's milk powder. *International Journal of Food Science and Nutrition Engineering*, 4, 15–19.